


Product Datasheet

TRIM38 antibody (orb1320922)

Description	Mouse monoclonal antibody to TRIM38
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-265 of human TRIM38 (NP_006346) produced in HEK293T cell.
Target	TRIM38
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53.2 kDa
Uniprot ID	O00635
NCBI	NM_006355
Dilution Range	WB 1:200~4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...