


Product Datasheet

LZIC antibody (orb1320990)

Description	Mouse monoclonal antibody to LZIC
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human LZIC(NP_115744) produced in HEK293T cell.
Target	LZIC
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.3 kDa
Uniprot ID	Q8WZA0
NCBI	NM_032368
Dilution Range	WB 1:200~4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...