

Product Datasheet

NRAS antibody (orb1321040)

Description	Mouse monoclonal antibody to NRAS
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human NRAS (NP_002515) produced in HEK293T cell.
Target	NRAS
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21 kDa
Uniprot ID	P01111
NCBI	NM_002524
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...