


Product Datasheet

LZIC antibody (orb1321043)

Description	Mouse monoclonal antibody to LZIC
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human LZIC(NP_115744) produced in E.coli.
Target	LZIC
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	21.3 kDa
Uniprot ID	Q8WZA0
NCBI	NM_032368
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...