


Product Datasheet

TBC1D13 antibody (orb1321071)

Description	Mouse monoclonal antibody to TBC1D13
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TBC1D13 (NP_060671) produced in HEK293T cell.
Target	TBC1D13
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46.4 kDa
Uniprot ID	Q9NVG8
NCBI	NM_018201
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...