


Product Datasheet

BBOX1 antibody (orb1321076)

Description	Mouse monoclonal antibody to BBOX1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human BBOX1(NP_003977) produced in HEK293T cell.
Target	BBOX1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	44.5 kDa
Uniprot ID	O75936
NCBI	NM_003986
Dilution Range	WB 1:4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...