


Product Datasheet

MEK1 (MAP2K1) antibody (orb1321118)

Description	Mouse monoclonal antibody to MEK1 (MAP2K1)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human MAP2K1(NP_002746) produced in HEK293T cell.
Target	MAP2K1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	43.3 kDa
Uniprot ID	Q02750
NCBI	NM_002755
Dilution Range	WB 1:400~4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...