


Product Datasheet

ATXN7L1 antibody (orb1321135)

Description	Mouse monoclonal antibody to ATXN7L1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ATXN7L1 (NP_689962) produced in HEK293T cell.
Target	ATXN7L1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	16 kDa
Uniprot ID	Q9ULK2
NCBI	NM_152749
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...