

Product Datasheet

C10orf7 (CDC123) antibody (orb1321151)

Description Mouse monoclonal antibody to C10orf7 (CDC123)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen Full length human recombinant protein of human CDC123(NP_006014) produced in HEK293T cell.

Target CDC123

Note For research use only

Isotype IgG2b

Clonality Monoclonal


MW 39 kDa

Uniprot ID [O75794](#)

NCBI [NM_006023](#)

Dilution Range WB 1:4000, IHC 1:150, IF 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...