


Product Datasheet

RAB24 antibody (orb1321154)

Description	Mouse monoclonal antibody to RAB24
Species/Host	Mouse
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human RAB24(NP_001026847) produced in HEK293T cell.
Target	RAB24
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	22.9 kDa
Uniprot ID	Q969Q5
NCBI	NM_001031677
Dilution Range	WB 1:200~4000, IF: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...