


Product Datasheet

CD272 (BTLA) antibody (orb1321157)

Description	Mouse monoclonal antibody to CD272 (BTLA)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human BTLA(NP_861445) produced in HEK293T cell.
Target	BTLA
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	32.6 kDa
Uniprot ID	Q7Z6A9
NCBI	NM_181780
Dilution Range	WB 1:4000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...