


Product Datasheet

ADAMTS8 antibody (orb1321160)

Description	Mouse monoclonal antibody to ADAMTS8
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human ADAMTS8(NP_008968) produced in HEK293T cell.
Target	ADAMTS8
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	96.3 kDa
Uniprot ID	Q9UP79
NCBI	NM_007037
Dilution Range	WB 1:1000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...