


Product Datasheet

DIXDC1 antibody (orb1321223)

Description	Mouse monoclonal antibody to DIXDC1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human DIXDC1(NP_219493) produced in HEK293T cell.
Target	DIXDC1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	53.8 kDa
Uniprot ID	Q155Q3
NCBI	NM_033425
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...