


Product Datasheet

RHOJ antibody (orb1321235)

Description	Mouse monoclonal antibody to RHOJ
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human RHOJ(NP_065714) produced in HEK293T cell.
Target	RHOJ
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	23.6 kDa
Uniprot ID	Q9H4E5
NCBI	NM_020663
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...