


Product Datasheet

PGD antibody (orb1321243)

Description	Mouse monoclonal antibody to PGD
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human PGD(NP_002622) produced in HEK293T cell.
Target	PGD
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53 kDa
Uniprot ID	P52209
NCBI	NM_002631
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...