


Product Datasheet

Tapasin Related Protein (TAPBPL) antibody (orb1321276)

Description	Mouse monoclonal antibody to Tapasin Related Protein (TAPBPL)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TAPBPL(NP_060479) produced in HEK293T cell.
Target	TAPBPL
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	50 kDa
Uniprot ID	Q9BX59
NCBI	NM_018009
Dilution Range	WB 1:2000, IHC 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...