


Product Datasheet

CATIP antibody (orb1321300)

Description	Mouse monoclonal antibody to CATIP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human C2orf62(NP_940961) produced in HEK239T cell.
Target	CATIP
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	43.7 kDa
Uniprot ID	Q7Z7H3
NCBI	NM_198559
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...