


Product Datasheet

NDOR1 antibody (orb1321334)

Description	Mouse monoclonal antibody to NDOR1
Species/Host	Mouse
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human NDOR1(NP_055249) produced in HEK293T cell.
Target	NDOR1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	66.6 kDa
Uniprot ID	Q9UHB4
NCBI	NM_014434
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...