


Product Datasheet

KIF25 antibody (orb1321363)

Description	Mouse monoclonal antibody to KIF25
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human KIF25(NP_005346) produced in HEK293T cell.
Target	KIF25
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.1 kDa
Uniprot ID	Q9UIL4
NCBI	NM_005355
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...