


Product Datasheet

HTF9C (TRMT2A) antibody (orb1321372)

Description	Mouse monoclonal antibody to HTF9C (TRMT2A)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human TRMT2A(NP_073564) produced in HEK293T cell.
Target	TRMT2A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	68.5 kDa
Uniprot ID	Q8IZ69
NCBI	NM_022727
Dilution Range	WB 1:400~4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...