


Product Datasheet

ARMC1 antibody (orb1321373)

Description	Mouse monoclonal antibody to ARMC1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human ARMC1(NP_060590) produced in HEK293T cell.
Target	ARMC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	31.1 kDa
Uniprot ID	Q9NVT9
NCBI	NM_018120
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...