


Product Datasheet

CD1C antibody (orb1321384)

Description	Mouse monoclonal antibody to CD1C
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human CD1C(NP_001756) produced in HEK293T cell.
Target	CD1C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.7 kDa
Uniprot ID	P29017
NCBI	NM_001765
Dilution Range	WB 1:400~4000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...