

Product Datasheet

TSC22D3 antibody (orb1321399)

Description	Mouse monoclonal antibody to TSC22D3
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TSC22D3(NP_004080) produced in HEK293T cell.
Target	TSC22D3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	14.6 kDa
Uniprot ID	Q99576
NCBI	NM_004089
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...