

Product Datasheet

MAGEA4 antibody (orb1321410)

Description	Mouse monoclonal antibody to MAGEA4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human MAGEA4(NP_001011550) produced in HEK293T cell.
Target	MAGEA4
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	34.7 kDa
Uniprot ID	P43358
NCBI	NM_001011550
Dilution Range	WB 1:4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...