


Product Datasheet

PPIL3 antibody (orb1321417)

Description	Mouse monoclonal antibody to PPIL3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human PPIL3(NP_572028) produced in E.coli.
Target	PPIL3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	18 kDa
Uniprot ID	Q9H2H8
NCBI	NM_131916
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...