


Product Datasheet

RAB37 antibody (orb1321447)

Description	Mouse monoclonal antibody to RAB37
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human RAB37(NP_783865) produced in HEK293T cell.
Target	RAB37
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	24 kDa
Uniprot ID	Q96AX2
NCBI	NM_175738
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...