

Product Datasheet

FATE1 antibody (orb1321476)

Description	Mouse monoclonal antibody to FATE1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human FATE1(NP_055814) produced in HEK293T cell.
Target	FATE1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	20.5 kDa
Uniprot ID	Q969F0
NCBI	NM_033085
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...