


Product Datasheet

LRCH4 antibody (orb1321490)

Description	Mouse monoclonal antibody to LRCH4
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human LRCH4 (NP_002310) produced in HEK293T cell.
Target	LRCH4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	73.3 kDa
Uniprot ID	O75427
NCBI	NM_002319
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...