


Product Datasheet

CGI 62 (ZC2HC1A) antibody (orb1321544)

Description	Mouse monoclonal antibody to CGI 62 (ZC2HC1A)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human FAM164A(NP_057094) produced in E.coli.
Target	ZC2HC1A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	34.9 kDa
Uniprot ID	Q96GY0
NCBI	NM_016010
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...