


Product Datasheet

ALX4 antibody (orb1321546)

Description	Mouse monoclonal antibody to ALX4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human ALX4(NP_068745) produced in HEK293T cell.
Target	ALX4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	44.1 kDa
Uniprot ID	Q9H161
NCBI	NM_021926
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...