


Product Datasheet

NDUFB10 antibody (orb1321548)

Description	Mouse monoclonal antibody to NDUFB10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human NDUFB10(NP_004539) produced in E.coli.
Target	NDUFB10
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	20.6 kDa
Uniprot ID	O96000
NCBI	NM_004548
Dilution Range	WB 1:1000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...