


Product Datasheet

UBE2S antibody (orb1321550)

Description	Mouse monoclonal antibody to UBE2S
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human UBE2S(NP_055316) produced in E.coli.
Target	UBE2S
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	23.7 kDa
Uniprot ID	Q16763
NCBI	NM_014501
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...