


Product Datasheet

p15 INK4b (CDKN2B) antibody (orb1321552)

Description	Mouse monoclonal antibody to p15 INK4b (CDKN2B)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human CDKN2B(NP_004927) produced in HEK293T cell.
Target	CDKN2B
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	14.5 kDa
Uniprot ID	P42772
NCBI	NM_004936
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...