


Product Datasheet

PLPBP antibody (orb1321554)

Description	Mouse monoclonal antibody to PLPBP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PROSC(NP_009129) produced in E.coli.
Target	PLPBP
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	30.2 kDa
Uniprot ID	O94903
NCBI	NM_007198
Dilution Range	WB 1:200~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...