


Product Datasheet

PHF7 antibody (orb1321561)

Description	Mouse monoclonal antibody to PHF7
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 155-381 of human PHF7(NP_057567) produced in E.coli.
Target	PHF7
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43.6 kDa
Uniprot ID	Q9BWX1
NCBI	NM_016483
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...