


Product Datasheet

SLAP2 (SLA2) antibody (orb1321573)

Description	Mouse monoclonal antibody to SLAP2 (SLA2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SLA2(NP_115590) produced in HEK293T cell.
Target	SLA2
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	28.4 kDa
Uniprot ID	Q9H6Q3
NCBI	NM_032214
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...