


Product Datasheet

RIOK2 antibody (orb1321575)

Description	Mouse monoclonal antibody to RIOK2
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 277-552 of human RIOK2(NP_060813) produced in E.coli.
Target	RIOK2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	63.1 kDa
Uniprot ID	Q9BVS4
NCBI	NM_018343
Dilution Range	WB 1:500~4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...