


Product Datasheet

RAB20 antibody (orb1321607)

Description	Mouse monoclonal antibody to RAB20
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human RAB20(NP_060287) produced in HEK293T cell.
Target	RAB20
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	26.1 kDa
Uniprot ID	Q9NX57
NCBI	NM_017817
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...