


Product Datasheet

CHN 1 (CHN1) antibody (orb1321627)

Description	Mouse monoclonal antibody to CHN 1 (CHN1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human CHN1(NP_001813) produced in HEK293T cell.
Target	CHN1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	53 kDa
Uniprot ID	P15882
NCBI	NM_001822
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...