


Product Datasheet

TSC22D1 antibody (orb1321635)

Description	Mouse monoclonal antibody to TSC22D1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human TSC22D1(NP_006013) produced in HEK293T cell.
Target	TSC22D1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	15.5 kDa
Uniprot ID	Q15714
NCBI	NM_006022
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...