


Product Datasheet

DREF (ZBED1) antibody (orb1321653)

Description	Mouse monoclonal antibody to DREF (ZBED1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZBED1(NP_004720) produced in HEK293T cell.
Target	ZBED1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	78 kDa
Uniprot ID	O96006
NCBI	NM_004729
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...