

Product Datasheet

STING (TMEM173) antibody (orb1321668)

Description Mouse monoclonal antibody to STING (TMEM173)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human TMEM173(NP_938023) produced in HEK293T cell.

Target STING1

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 42 kDa

Uniprot ID [Q86WV6](#)

NCBI [NM_198282](#)

Dilution Range WB 1:200 - 1:1000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...