


Product Datasheet

MAFB antibody (orb1321675)

Description	Mouse monoclonal antibody to MAFB
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 32-320 of human MAFB (NP_005452) produced in E.coli.
Target	MAFB
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.6 kDa
Uniprot ID	Q9Y5Q3
NCBI	NM_005461
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...