


Product Datasheet

MMACHC antibody (orb1321685)

Description	Mouse monoclonal antibody to MMACHC
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human MMACHC(NP_056321) produced in HEK293T cell.
Target	MMACHC
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	31.5 kDa
Uniprot ID	Q9Y4U1
NCBI	NM_015506
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...