


Product Datasheet

Tumor protein D52 like 3 (TPD52L3) antibody (orb1321705)

Description	Mouse monoclonal antibody to Tumor protein D52 like 3 (TPD52L3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TPD52L3(NP_001001875) produced in HEK293T cell.
Target	TPD52L3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	14.4 kDa
Uniprot ID	Q96J77
NCBI	NM_001001875
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...