


Product Datasheet

C20orf30 (TMEM230) antibody (orb1321723)

Description	Mouse monoclonal antibody to C20orf30 (TMEM230)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human C20orf30(NP_054864) produced in HEK293T cell.
Target	TMEM230
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	13 kDa
Uniprot ID	Q96A57
NCBI	NM_014145
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...