


Product Datasheet

COMMD1 antibody (orb1321753)

Description	Mouse monoclonal antibody to COMMD1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human COMMD1(NP_689729) produced in HEK293T cell.
Target	COMMD1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	21 kDa
Uniprot ID	Q8N668
NCBI	NM_152516
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...