


Product Datasheet

AFG3L2 antibody (orb1321787)

Description	Mouse monoclonal antibody to AFG3L2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protien fragment corresponding to amino acids 1-250 of human AFG3L2(NP_006787) produced in E.coli.
Target	AFG3L2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	88.3 kDa
Uniprot ID	Q9Y4W6
NCBI	NM_006796
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...