


Product Datasheet

SENP2 antibody (orb1321793)

Description	Mouse monoclonal antibody to SENP2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 139-523 of human SENP2(NP_067640) produced in E.coli.
Target	SENP2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q9HC62
NCBI	NM_021627
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...