


Product Datasheet

SENP2 antibody (orb1321817)

Description	Mouse monoclonal antibody to SENP2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 139-523 of human SENP2(NP_067640) produced in E.coli.
Target	SENP2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9HC62
NCBI	NM_021627
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot
analysis of
HEK293T
cells w...